




Cat No: HRC042 – HER2 IHC control Array

Lot#	Cores	Size	Cut	Format	QA/QC
HRC04201	4	4.5mm	4um	1X4	H&E, IHC anti-HER2



Recommended applications: HER2 IHC control with strong, moderate, low and negative expressers

Description: The array contains four cases of invasive ductal carcinoma of the breast. All the tissues were from surgical resection. They were fixed in 10% neutral buffered formalin for 24 hours and processed using identical SOPs. Sections were picked onto Superfrost Plus or Startfrost slides. As validated by IHC, A1 is a strong (+++) expresser, A2 is a moderate (++) expresser, A3 is a weak expresser for rabbit monoclonal HER2 antibody, while A4 is a weak expresser for DAKO's HER2 antibody.

Notes: Bake at 60C for 30 minutes before use. If antigen retrieving is needed, it is always a good idea to start with a protocol with weak to mild strength

Certified by: Langxing Pan, M.D., Ph.D.