



Cat No: **LYM1501 - Lymphoma tissue array**

Lot#	Cores	Size	Cut	Format	QA/QC		
LYM150101	150	1.1mm	4um	10X15	H&E, IHC anti-LCA		

Recommended applications: For research use only. RNA or protein analyses using IHC or ISH suitable for general lymphoma tissue profiling.

Description: Lymphoma array, 75 cases of various types of nodal and extranodal lymphoma in duplicates. All the tissues were from surgical resection. They were fixed in 10% neutral buffered formalin for 24 hours and processed using identical SOPs. Sections were picked onto Superfrost Plus or APES coated Superfrost slides. They all have a guaranteed six months' shelf-life at 4C from the date of shipment. Each slide has >95% tissue core retention.

Array position	Age	Sex	Anatomic site	Histology
A1,B1	26	M	Lymph node	Reactive
A2,B2	47	M	Lymph node	Reactive
A3,B3	31	M	Tonsil	Reactive
A4,B4	14	M	Tonsil	Reactive
A5,B5	21	F	Spleen	Normal
A6,B6	77	F	Thyroid	MALT B cell lymphoma
A7,B7	62	M	Stomach	Non-Hodgkin B cell lymphoma
A8,B8	53	M	Small intestine	Peripheral T cell lymphoma, unspecified
A9,B9	41	M	Neck lymph node	Hodgkin lymphoma
A10,B10	55	M	Neck lymph node	Peripheral T cell lymphoma, unspecified?
A11,B11	54	M	Shoulder	Non-Hodgkin B cell lymphoma
A12,B12	29	M	Neck lymph node	Non-Hodgkin B cell lymphoma
A13,B13	5	M	Lymph node	Anaplastic large cell lymphoma?
A14,B14	40	M	Lymph node	Hodgkin lymphoma
A15,B15	48	M	Colon	Non-Hodgkin B cell lymphoma
C1,D1	60	M	Colon	Non-Hodgkin B cell lymphoma
C2,D2	42	M	Axillary lymph node	Non-Hodgkin B cell lymphoma
C3,D3	7	M	Lymph node	Hodgkin lymphoma
C4,D4	66	M	Lymph node	Non-Hodgkin lymphoma
C5,D5	40	M	Ilium-cecum	Non-Hodgkin B cell lymphoma
C6,D6	33	M	Neck lymph node	Non-Hodgkin B cell lymphoma
C7,D7	51	M	Neck lymph node	Anaplastic large cell lymphoma?
C8,D8	48	M	Lung lymph node	Non-Hodgkin B cell lymphoma?
C9,D9	72	M	Testis	Non-Hodgkin B cell lymphoma
C10,D10	71	F	Lymph node	Non-Hodgkin B cell lymphoma
C11,D11	43	M	Stomach	MALT B cell lymphoma, high grade
C12,D12	52	F	Neck lymph node	Diffuse large B cell lymphoma
C13,D13	60	M	Colon	MALT B cell lymphoma?
C14,D14	30	M	Peritoneal cavity	Hodgkin lymphoma
C15,D15	40	F	Peritoneal cavity	Non-Hodgkin B cell lymphoma
E1,F1	33	F	Stomach	MALT B cell lymphoma, high grade
E2,F2	70	F	Mesenteric lymph node	Non-Hodgkin B cell lymphoma
E3,F3	43	F	Neck lymph node	Non-Hodgkin B cell lymphoma
E4,F4	65	M	Lymph node	Non-Hodgkin B cell lymphoma
E5,F5	76	M	Tonsil	Diffuse large B cell lymphoma
E6,F6	32	F	Pelvic Cavity	Non-Hodgkin B cell lymphoma



Advancing Biomedical Science Through Tissue Arrays

E7,F7	39	M	Neck lymph node	Non-Hodgkin B cell lymphoma
E8,F8	50	M	Neck lymph node	Hodgkin lymphoma
E9,F9	43	M	Ilium	Non-Hodgkin lymphoma
E10,F10	53	F	Neck lymph node	Mantle cell lymphoma
E11,F11	72	F	Lymph node	Non-Hodgkin B cell lymphoma
E12,F12	53	M	Neck lymph node	Diffuse large B cell lymphoma
E13,F13	55	F	Mesenteric lymph node	Non-Hodgkin lymphoma
E14,F14	55	F	Breast	Non-Hodgkin B cell lymphoma
E15,F15	13	M	Shoulder	Diffuse large B cell lymphoma
G1,H1	35	M	Lymph node	Non-Hodgkin B cell lymphoma
G2,H2	40	M	Tonsil	Non-Hodgkin B cell lymphoma
G3,H3	54	F	Lymph node	Hodgkin lymphoma
G4,H4	24	M	Ilium-cecum	Non-Hodgkin B cell lymphoma
G5,H5	64	M	Lymph node	Non-Hodgkin B cell lymphoma
G6,H6	49	M	Neck lymph node	Hodgkin lymphoma
G7,H7	65	M	Stomach	MALT B cell lymphoma, high grade
G8,H8	73	M	Kidney	Non-Hodgkin B cell lymphoma
G9,H9	4	M	Neck lymph node	Non-Hodgkin B cell lymphoma
G10,H10	70	F	Ilium	Non-Hodgkin B cell lymphoma
G11,H11	52	F	Lung	Non-Hodgkin lymphoma
G12,H12	51	M	Stomach	Non-Hodgkin B cell lymphoma
G13,H13	28	M	Ilium	Non-Hodgkin B cell lymphoma
G14,H14	58	F	Stomach	MALT B cell lymphoma, high grade
G15,H15	38	M	Lymph node	Non-Hodgkin B cell lymphoma
I1,J1	58	F	Stomach	MALT B cell lymphoma, high grade
I2,J2	67	M	Lymph node	Non-Hodgkin lymphoma
I3,J3	53	F	Stomach	MALT B cell lymphoma, high grade
I4,J4	73	M	Lymph node	Hodgkin lymphoma
I5,J5	17	M	Small intestine	Non-Hodgkin B cell lymphoma
I6,J6	43	M	Brain	Non-Hodgkin B cell lymphoma
I7,J7	43	M	Neck lymph node	Non-Hodgkin B cell lymphoma
I8,J8	76	M	Neck lymph node	Non-Hodgkin B cell lymphoma
I9,J9	51	F	Neck lymph node	Hodgkin lymphoma
I10,J10	74	F	Neck lymph node	Non-Hodgkin lymphoma
I11,J11	57	M	Neck lymph node	Hodgkin lymphoma
I12,J12	42	F	Breast	Non-Hodgkin B cell lymphoma
I13,J13	42	F	Neck lymph node	Hodgkin lymphoma
I14,J14	38	F	Neck lymph node	Hodgkin lymphoma
I15,J15	44	F	Peritoneal cavity	Anaplastic large cell lymphoma

Notes: Bake at 60C for 30 minutes before use. If antigen retrieving is needed, it is always a good idea to start with a protocol with weak to mild strength.

Certified by: Fancai Li, M.D.