Phone: 415-863-2380

Support: <u>info@pantomics.com</u> Web: www.pantomics.com

Advancing Biomedical Science Through Tissue Analysis

# Cat No: NPC1261 - Nasopharyngeal carcinoma tissue array

Lot#	Cores	Size	Cut	Format	QA/QC	J J	1 2 3 4 5 6 7 8 9 10 11 12 13 14
NPC126101	126	1.1mm	4um		H&E, IHC anti- Cytokeratin	Pantomics, In	B B B B B B B B B B B B B B B B B B B

**Recommended Application:** For research use only. Designed for IHC or ISH based protein or RNA tissue profiling in nasopharyngeal carcinoma.

**Description:** Nasopharyngeal carcinoma tissue array, set 1, non-overlapping with other NPC tissue arrays, containing 96 cases of NPC and 30 cases of normal-reactive nasopharyngeal mucosal tissues.

All the tissues were from surgical resection. They were fixed in 10% neutral buffered formalin for 24 hours and processed using identical SOPs. Sections were picked onto Superfrost Plus or APES coated Superfrost slides. They can be stored for use at 4C for up to six months from the date of shipment. There may be 5 to 10% of tissue core loss.

Array Position	Sex	Age	Anatomic Site	Pathology	Grade	Stage
A01	М	62	Head and neck, nasopharynx	Chronic inflammation	null	null
A02	М	50	Head and neck, nasopharynx	Chronic inflammation	null	null
A03	М	49	Head and neck, nasopharynx	Chronic inflammation	null	null
A04	М	42	Head and neck, nasopharynx	Chronic inflammation	null	null
A05	М	37	Head and neck, nasopharynx	Chronic inflammation	null	null
A06	М	46	Head and neck, nasopharynx	Chronic inflammation	null	null
A07	М	45	Head and neck, nasopharynx	Chronic inflammation	null	null
A08	М	59	Head and neck, nasopharynx	Chronic inflammation	null	null
A09	М	72	Head and neck, nasopharynx	Chronic inflammation	null	null
A10	F	23	Head and neck, nasopharynx	Chronic inflammation	null	null
A11	F	59	Head and neck, nasopharynx	Chronic inflammation	null	null
A12	М	36	Head and neck, nasopharynx	Chronic inflammation	null	null
A13	М	18	Head and neck, nasopharynx	Chronic inflammation	null	null
A14	F	54	Head and neck, nasopharynx	Chronic inflammation	null	null
B01	F	13	Head and neck, nasopharynx	Chronic inflammation	null	null
B02	F	49	Head and neck, nasopharynx	Chronic inflammation	null	null
B03	F	39	Head and neck, nasopharynx	Chronic inflammation	null	null
B04	М	29	Head and neck, nasopharynx	Chronic inflammation	null	null
B05	F	63	Head and neck, nasopharynx	Chronic inflammation	null	null
B06	М	50	Head and neck, nasopharynx	Chronic inflammation	null	null
B07	F	58	Head and neck, nasopharynx	Chronic inflammation	null	null
B08	F	58	Head and neck, nasopharynx	Chronic inflammation	null	null
B09	М	43	Head and neck, nasopharynx	Chronic inflammation	null	null
B10	М	55	Head and neck, nasopharynx	Chronic inflammation	null	null
B11	М	20	Head and neck, nasopharynx	Chronic inflammation	null	null
B12	М	42	Head and neck, nasopharynx	Chronic inflammation	null	null
B13	F	38	Head and neck, nasopharynx	Chronic inflammation	null	null
B14	F	37	Head and neck, nasopharynx	Chronic inflammation	null	null
C01	F	50	Head and neck, nasopharynx	Chronic inflammation	null	null
C02	М	27	Head and neck, nasopharynx	Chronic inflammation	null	null
C03	М	33	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C04	F	37	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C05	М	45	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C06	М	36	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C07	М	63	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null

Pantomics, Inc., 1049 Union Ave, Suite D, CA 94533 Phone: 415-863-2380; Fax: 510-225-1084; Email: info@pantomics.com

Phone: 415-863-2380 Support: <u>info@pantomics.com</u> Web: www.pantomics.com

			ics/ Inc	Web. www.pe		
		Adv	ancing Biomedical Science	e Through Tissue Amaya		
C08	М	48	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C09	М	48	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C10	F	43	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C11	М	19	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C12	М	49	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C13	М	41	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
C14	М	43	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D01	F	58	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D02	F	44	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D03	М	67	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D04	М	25	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D05	М	62	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D06	F	41	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D07	F	50	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D08	М	60	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D09	F	27	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D10	М	48	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D11	М	46	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D12	М	59	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D13	М	60	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
D14	F	44	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
E01	F	67	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
E02	М	65	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
E03	М	45	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
E04	М	50	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
E05	М	49	Head and neck, nasopharynx	Non-keratinizing undifferentiated	III	null
E06	М	44	Head and neck, nasopharynx	Carcinoma  Non-keratinizing undifferentiated	III	null
E07	М	45	Head and neck, nasopharynx	Carcinoma  Non-keratinizing undifferentiated	III	null
E08	М	37	Head and neck, nasopharynx	carcinoma  Non-keratinizing undifferentiated	III	null
E09	М	52	Head and neck, nasopharynx	Carcinoma  Non-keratinizing undifferentiated	III	null
E10	F	40	Head and neck, nasopharynx	Carcinoma  Non-keratinizing undifferentiated	III	null
E11	М	47	Head and neck, nasopharynx	carcinoma  Non-keratinizing undifferentiated	III	null
E12	М	46	Head and neck, nasopharynx	carcinoma  Non-keratinizing undifferentiated	III	null
CIZ				carcinoma		

Phone: 415-863-2380

Support: <u>info@pantomics.com</u> Web: www.pantomics.com

		Adv	ancing Biomedical Scienc	e Through Tissue Amana		
************	*****			carcinoma		
E14	М	46	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F01	М	62	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F02	М	50	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F03	М	55	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F04	М	57	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F05	М	59	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F06	F	64	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F07	М	59	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F08	F	41	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F09	М	50	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F10	F	40	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F11	М	32	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F12	М	54	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F13	F	64	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
F14	М	59	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G01	М	60	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G02	М	50	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G03	М	53	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G04	М	49	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G05	М	49	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G06	М	62	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G07	М	60	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G08	М	42	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G09	М	57	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G10	М	57	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G11	М	39	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G12	F	50	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G13	М	60	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
G14	М	48	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H01	М	42	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H02	F	50	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H03	М	61	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H04	М	54	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null

Phone: 415-863-2380 Support: <u>info@pantomics.com</u>

Web: www.pantomics.com

000000000000000000000000000000000000000		STREET, 477	THE RESERVE OF THE PARTY OF THE	W-1	THE REAL PROPERTY.	
		Adv	ancing Biomedical Scienc			
H05	М	47	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H06	М	49	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H07	М	56	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H08	F	53	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H09	М	51	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H10	М	36	Head and neck, nasopharynx	Non-keratinizing undifferentiated carcinoma	III	null
H11	F	48	Head and neck, nasopharynx	Non-keratinizing undifferentiated	III	null
H12	М	48	Head and neck, nasopharynx	carcinoma   Non-keratinizing undifferentiated	III	null
H13	М	59	Head and neck, nasopharynx	carcinoma   Non-keratinizing undifferentiated	III	null
H14	М	46	Head and neck, nasopharynx	carcinoma   Non-keratinizing undifferentiated	III	null
I01	F	58	Head and neck, nasopharynx	Non-keratinizing undifferentiated	III	null
I02	М	19	Head and neck, nasopharynx	Non-keratinizing undifferentiated	III	null
I03	М	77	Head and neck, nasopharynx	carcinoma   Non-keratinizing undifferentiated	III	null
I04	М	43	Head and neck, nasopharynx	carcinoma   Non-keratinizing undifferentiated	III	null
I05	F	33	Head and neck, nasopharynx	carcinoma   Non-keratinizing undifferentiated	III	null
106	F	44	Head and neck, nasopharynx	Non-keratinizing undifferentiated	III	null
107	F	48	Head and neck, nasopharynx	Non-keratinizing undifferentiated	III	null
108	М	52	Head and neck, nasopharynx	Non-keratinizing undifferentiated	III	null
109	М	53	Head and neck, nasopharynx	carcinoma   Non-keratinizing undifferentiated	III	null
I10	F	57	Head and neck, nasopharynx	carcinoma Non-keratinizing undifferentiated	III	null
I11	М	36	Head and neck, nasopharynx	Non-keratinizing undifferentiated	III	null
I12	М	48	Head and neck, nasopharynx	carcinoma   Non-keratinizing undifferentiated	III	null
I13	М	29	Head and neck, nasopharynx	Non-keratinizing undifferentiated	III	null
I14	M	64	Head and neck, nasopharynx	carcinoma   Non-keratinizing undifferentiated   carcinoma	III	null

**Notes:** Bake at 60C for ~60 minutes before use. If antigen retrieving is needed, it is important to avoid **direct-boiling and high pH or high strength** antigen retrieving buffer. For availability of complimentary IHC data, please contact us at info@pantomics.com.

Certified by: Fancai Li, M.D.

Phone: 415-863-2380

Support: <u>info@pantomics.com</u>
Web: www.pantomics.com

Advancing Biomedical Science Through Tissue Ama

# TNM Classification: Nasopharyngeal carcinoma

### **T - Primary Tumor**

- TX Primary tumor cannot be assessed
- T0 No evidence of primary tumor
- Tis Carcinoma in situ
- T1 Tumor confined to nasopharynx
- T2 Tumor extends to soft tissue of oropharynx and/or nasal fossa
- T2a without parapharyngeal extension
- T2b with parapharyngeal extension
- T3 Tumor invades bony structures and/or paranasal sinuses
- T4 Tumor with intracranial extension and/or involvement of cranial nerves, infratemporal fossa, hypopharynx, or orbit

### N - Regional Lymph nodes

- NX Reginal lymph nodes cannot be assessed
- NO No regional lymph node metastasis
- N1 Unilateral metastasis in lymph node(s), 6cm or less in greatest dimension, above supraclavicular fossa
- N2 Bilateral metastasis in lymph node(s), 6cm or less in greatest dimension, above supraclavicular fossa
- N3 Metastasis in lymph node(s)
  - (a) greater than 6cm in dimension
  - (b) in the supraclavicular

#### **M** - Distant Metastasis

- MX Distant metastasis cannot be assessed
- M0 No distant metastasis
- M1 Distant metastasis