Phone: 415-863-2380

Support: <u>info@pantomics.com</u>
Web: www.pantomics.com

Advancing Biomedical Science Through Tissue Anaya:

Cat No: **LYM1504**- **Lymphoma tissue array**

Lot#	Cores	Size	Cut	Format	QA/QC	1 2 3 4 5 6 7 8 9 10 11 12 13 14	15
LYM150401	150	1.1mm	4um	10X15	H&E, IHC anti-LCA	Pantomics, Inc	+

Recommended applications: For research use only. For RNA or protein profiling of general lymphoma/Non-Hodgkin B cell lymphoma tissue using IHC or ISH, and for antibody characterization.

Description: Lymphoma array, 75 cases, set 4, non-overlapping with LYM1501, LYM1502, LYM1503 or LYM1505, 150 cores from reactive lymphoid tissues (5 cases) and various types of nodal and extranodal lymphoma (70 cases) in duplicates.

All the tissues were from surgical resection. They were fixed in 10% neutral buffered formalin for 24 hours and processed using identical SOPs. Sections were picked onto Superfrost Plus or APES coated Superfrost slides. They can be stored for use at 4C for up to six months from the date of shipment. **There may be 5 to 10% of tissue core loss**.

Array Position	Sex	Age	Anatomic Site	Pathology	Grade	Stage
A01	М	40	Tonsil	Inflammation, chronic		
A02	М	31	Tonsil	Inflammation, chronic		
A03	М	40	Spleen	Congestion		
A04	F	69	Lymph node	Normal		
A05	F	46	Lymph node	Normal		
A06	F	69	Tonsil	Lymphoma, diffuse large B-cell lymphoma		
A07	М	35	Lymph node	Lymphoma, lymph plasmacytoid lymphoma		
A08	М	42	Lymph node, inguinal	Lymphoma, diffuse large B-cell lymphoma		
A09	F	33	Mediastinum	Lymphoma, diffuse large B-cell lymphoma		
A10	М	47	Lymph node, neck	Lymphoma, non-Hodgkin B-cell lymphoma		
A11	F	25	Tonsil	Lymphoma, non-Hodgkin B-cell lymphoma		
A12	F	17	Mesenteric site	Lymphoma, diffuse large B-cell lymphoma		
A13	М	66	Stomach	Lymphoma, non-Hodgkin B-cell lymphoma		
A14	F	63	Lymph node, neck	Lymphoma, Hodgkin lymphoma		
A15	F	64	Lymph node, inguinal	Lymphoma, non-Hodgkin lymphoma		
B01	М	40	Tonsil	Inflammation, chronic		
B02	М	31	Tonsil	Inflammation, chronic		
B03	М	40	Spleen	Congestion		
B04	F	69	Lymph node	Normal		
B05	F	46	Lymph node	Normal		
B06	F	69	Tonsil	Lymphoma, diffuse large B-cell lymphoma		
B07	М	35	Lymph node	Lymphoma, lymph plasmacytoid lymphoma		
B08	М	42	Lymph node, inguinal	Lymphoma, diffuse large B-cell lymphoma		
B09	F	33	Mediastinum	Lymphoma, diffuse large B-cell lymphoma		
B10	М	47	Lymph node, neck	Lymphoma, non-Hodgkin B-cell lymphoma		
B11	F	25	Tonsil	Lymphoma, non-Hodgkin B-cell lymphoma		
B12	F	17	Mesenteric site	Lymphoma, diffuse large B-cell lymphoma		
B13	М	66	Stomach	Lymphoma, non-Hodgkin B-cell lymphoma		
B14	F	63	Lymph node, neck	Lymphoma, Hodgkin lymphoma		
B15	F	64	Lymph node, inguinal	Lymphoma, non-Hodgkin lymphoma		
C01	F	54	Lymph node	Lymphoma, follicular lymphoma		
C02	F	70	Head and neck	Lymphoma, diffuse large B-cell lymphoma		
C03	М	75	Head and neck, salivary gland	Lymphoma, non-Hodgkin B-cell lymphoma		
C04	М	67	Lymph node	Lymphoma, diffuse large B-cell lymphoma		

Pantomics, Inc., 1049 Union Ave, Suite D, CA 94533 Phone: 415-863-2380; Fax: 510-225-1084; Email: info@pantomics.com

Phone: 415-863-2380

Support: <u>info@pantomics.com</u> Web: www.pantomics.com

000000000000000000000000000000000000000			ics/ Inc	Web. www.pantonnes.com	
COF	1 14			Science Through Tissue Amaye	
C05	M	74	Lymph node, inguinal	Lymphoma, B small lymphocytic lymphoma	
C06	M	71	Lymph node, inguinal	Lymphoma, non-Hodgkin lymphoma	
C07	M	59	Abdomen wall	Lymphoma, diffuse large B-cell lymphoma	
C08	M F	49	Stomach	Lymphoma, non-Hodgkin lymphoma	
C09		71	Lymph node	Lymphoma, diffuse large B-cell lymphoma	
C10	M	73	Eye	Lymphoma, MALT B cell lymphoma	
C11	F	55	Ovary	Lymphoma, diffuse large B-cell lymphoma	
C12	M	60	Head and neck	Lymphoma, diffuse large B-cell lymphoma	
C13	F	65	Spleen	Lymphoma, follicular lymphoma	
C14	F	61	Lymph node, inguinal	Lymphoma, non-Hodgkin T-cell lymphoma, NOS	
C15	M	59	Stomach	Lymphoma, non-Hodgkin B-cell lymphoma	
D01	F	54	Lymph node	Lymphoma, follicular lymphoma	
D02	F	70	Head and neck	Lymphoma, diffuse large B-cell lymphoma	
D03	М	75	Head and neck, salivary gland	Lymphoma, non-Hodgkin B-cell lymphoma	
D04	М	67	Lymph node	Lymphoma, diffuse large B-cell lymphoma	
D05	М	74	Lymph node, inguinal	Lymphoma, B small lymphocytic lymphoma	
D06	М	71	Lymph node, inguinal	Lymphoma, non-Hodgkin lymphoma	
D07	М	59	Abdomen wall	Lymphoma, diffuse large B-cell lymphoma	
D08	М	49	Stomach	Lymphoma, non-Hodgkin lymphoma	
D09	F	71	Lymph node	Lymphoma, diffuse large B-cell lymphoma	
D10	М	73	Eye	Lymphoma, MALT B cell lymphoma	
D11	F	55	Ovary	Lymphoma, diffuse large B-cell lymphoma	
D12	М	60	Head and neck	Lymphoma, diffuse large B-cell lymphoma	
D13	F	65	Spleen	Lymphoma, follicular lymphoma	
D14	F	61	Lymph node, inguinal	Lymphoma, non-Hodgkin T-cell lymphoma, NOS	
D15	М	59	Stomach	Lymphoma, non-Hodgkin B-cell lymphoma	
E01	М	45	Lymph node	Lymphoma, follicular lymphoma	
E02	М	43	Chest wall	non-Hodgkin lymphoma?	
E03	F	67	Lymph node, axillary	Lymphoma, diffuse large B-cell lymphoma	
E04	F	65	Lymph node, axillary	Lymphoma, non-Hodgkin T-cell lymphoma, NOS	
E05	F	20	Head and neck	Lymphoma, Hodgkin lymphoma	
E06	М	60	Lymph node	Lymphoma, non-Hodgkin B-cell lymphoma	
E07	М	70	Lymph node, neck	Lymphoma, B small lymphocytic lymphoma	
E08	М	84	Spleen	Lymphoma, non-Hodgkin B-cell lymphoma	
E09	М	23	Tonsil	Lymphoma, diffuse large B-cell lymphoma	
E10	М	4	Lymph node, neck	Lymphoma, Hodgkin lymphoma	
E11	М	85	Lymph node, neck	Lymphoma, non-Hodgkin B-cell lymphoma	
E12	F	52	Head and neck	Lymphoma, Hodgkin lymphoma	
E13	F	70	Head and neck	Lymphoma, mantle cell lymphoma	
E14	М	10	Lymph node, neck	Lymphoma, anaplastic lymphoma, ALK+	
E15	М	15	Head and neck	Lymphoma, diffuse large B-cell lymphoma	
F01	М	45	Lymph node	Lymphoma, follicular lymphoma	
F02	М	43	Chest wall	non-Hodgkin lymphoma?	
F03	F	67	Lymph node, axillary	Lymphoma, diffuse large B-cell lymphoma	
F04	F	65	Lymph node, axillary	Lymphoma, non-Hodgkin T-cell lymphoma, NOS	
F05	F	20	Head and neck	Lymphoma, Hodgkin lymphoma	
F06	М	60	Lymph node	Lymphoma, non-Hodgkin B-cell lymphoma	
F07	М	70	Lymph node, neck	Lymphoma, B small lymphocytic lymphoma	
F08	М	84	Spleen	Lymphoma, non-Hodgkin B-cell lymphoma	
F09	М	23	Tonsil	Lymphoma, diffuse large B-cell lymphoma	
F10	М	4	Lymph node, neck	Lymphoma, Hodgkin lymphoma	
F11	М	85	Lymph node, neck	Lymphoma, non-Hodgkin B-cell lymphoma	
		52	Head and neck	Lymphoma, Hodgkin lymphoma	

Phone: 415-863-2380

Support: <u>info@pantomics.com</u> Web: www.pantomics.com

			ies, iiie	Web. www.pantonnes.com	
		Adv		Science Through Tissue Amayu	
F13	F	70	Head and neck	Lymphoma, mantle cell lymphoma	
F14	М	10	Lymph node, neck	Lymphoma, anaplastic lymphoma, ALK+	
F15	М	15	Head and neck	Lymphoma, diffuse large B-cell lymphoma	
G01	М	62	Head and neck	Lymphoma, diffuse large B-cell lymphoma	
G02	М	60	Leg	Lymphoma, diffuse large B-cell lymphoma	
G03	М	57	Brain	Lymphoma, diffuse large B-cell lymphoma	
G04	М	60	Lymph node	Lymphoma, Angioimmunoblastic T-cell	
				lymphoma?	
G05	М	38	Head and neck	Lymphoma, follicular lymphoma	
G06	М	65	Lymph node	Lymphoma, diffuse large B-cell lymphoma	
G07	М	32	Lymph node, neck	Lymphoma, non-Hodgkin lymphoma	
G08	М	56	Lymph node	Lymphoma, Hodgkin lymphoma	
G09	М	80	Trunk	Lymphoma, non-Hodgkin lymphoma	
G10	М	80	Lymph node, axillary	Lymphoma, non-Hodgkin lymphoma	
G11	М	42	Tonsil	Lymphoma, follicular lymphoma	
G12	F	58	Stomach	Lymphoma, diffuse large B-cell lymphoma	
G13	F	47	Lung	Lymphoma, MALT B cell lymphoma	
G14	F	35	Mediastinum	Lymphoma, Hodgkin lymphoma	
G15	F	57	Stomach	Lymphoma, diffuse large B-cell lymphoma	
H01	M	62	Head and neck	Lymphoma, diffuse large B-cell lymphoma	
H02	М	60	Leg	Lymphoma, diffuse large B-cell lymphoma	
H03	М	57	Brain	Lymphoma, diffuse large B-cell lymphoma	
H04	М	60	Lymph node	Lymphoma, Angioimmunoblastic T-cell lymphoma?	
H05	М	38	Head and neck	Lymphoma, follicular lymphoma	
H06	М	65	Lymph node	Lymphoma, diffuse large B-cell lymphoma	
H07	М	32	Lymph node, neck	Lymphoma, non-Hodgkin lymphoma	
H08	М	56	Lymph node	Lymphoma, Hodgkin lymphoma	
H09	М	80	Trunk	Lymphoma, non-Hodgkin lymphoma	
H10	М	80	Lymph node, axillary	Lymphoma, non-Hodgkin lymphoma	
H11	М	42	Tonsil	Lymphoma, follicular lymphoma	
H12	F	58	Stomach	Lymphoma, diffuse large B-cell lymphoma	
H13	F	47	Lung	Lymphoma, MALT B cell lymphoma	
H14	F	35	Mediastinum	Lymphoma, Hodgkin lymphoma	
H15	F	57	Stomach	Lymphoma, diffuse large B-cell lymphoma	
I01	F	74	Skin	Lymphoma, diffuse large B-cell lymphoma	
101	F	61	Trunk	Lymphoma, non-Hodgkin lymphoma	
		79		Lymphoma, B small lymphocytic lymphoma	
I03	М		Lymph node, neck		
104	F	54	Lymph node	Lymphoma, B small lymphocytic lymphoma	
105	M	16	Lymph node	Lymphoma, Hodgkin lymphoma	
106	F	38	Thyroid	Lymphoma, MALT B cell lymphoma	
107	М	45	Lymph node	Lymphoma, mantle cell lymphoma	
108	М	48	Head and neck	Lymphoma, diffuse large B-cell lymphoma	
109	М	81	Lymph node, inguinal	Lymphoma, diffuse large B-cell lymphoma	
I10	М	61	Lymph node	Lymphoma, follicular lymphoma	
I11	F	58	Head and neck	Lymphoma, non-Hodgkin lymphoma	
I12	М	68	Lymph node, neck	Lymphoma, Hodgkin lymphoma	
I13	F	61	Lymph node, neck	Lymphoma, non-Hodgkin lymphoma	
I14	М	70	Intestine, colon	Lymphoma, diffuse large B-cell lymphoma	
I15	М	47	Spleen	Lymphoma, follicular lymphoma	
J01	F	74	Skin	Lymphoma, diffuse large B-cell lymphoma	
J02	F	61	Trunk	Lymphoma, non-Hodgkin lymphoma	
J03	M	79	Lymph node, neck	Lymphoma, B small lymphocytic lymphoma	
J04	F	54	Lymph node	Lymphoma, B small lymphocytic lymphoma	
J05	M	16	Lymph node	Lymphoma, Hodgkin lymphoma	
105	IVI	10	Lymph node	шунірноніа, пои <u>чукін і</u> унірноніа	

Phone: 415-863-2380 Support: <u>info@pantomics.com</u>

Web: www.pantomics.com

		Adv	vancing Biomedical S	cience Through Tissue Amays	
J06	F	38	Thyroid	Lymphoma, MALT B cell lymphoma	
J07	М	45	Lymph node	Lymphoma, mantle cell lymphoma	
J08	М	48	Head and neck	Lymphoma, diffuse large B-cell lymphoma	
J09	М	81	Lymph node, inguinal	Lymphoma, diffuse large B-cell lymphoma	
J10	М	61	Lymph node	Lymphoma, follicular lymphoma	
J11	F	58	Head and neck	Lymphoma, non-Hodgkin lymphoma	
J12	М	68	Lymph node, neck	Lymphoma, Hodgkin lymphoma	
J13	F	61	Lymph node, neck	Lymphoma, non-Hodgkin lymphoma	
J14	М	70	Intestine, colon	Lymphoma, diffuse large B-cell lymphoma	
J15	М	47	Spleen	Lymphoma, follicular lymphoma	

Notes: Bake at 60C for ~60 minutes before use. If antigen retrieving is needed, it is important to avoid **direct-boiling and high pH or high strength** antigen retrieving buffer. For availability of complimentary IHC data, please contact us at info@pantomics.com.

Certified by: Fancai Li, M.D.

Clinical Staging: Malignant lymphoma

- Stage I Involvement of a single lymph node region (I), or localized involvement of a single extrahepatic organ or site (I_E)
- Stage II Involvement of two or more lymph node regions on the same site of the diaphragm (II), or localized involvement of a single extralymphatic organ or site and its regional lymph node(s) with or without involvement of other lymph node regions on the same side of the diaphragm (II_E). The number of lymph node regions involved may be indicated by a subscript (e.g. II_3)
- Stage III Involvement of lymph node regions on both sides of the diaphragm (III), which may also be accompanied by localized involvement of an associated extralymphatic organ or sites (III_E), or by involvement of spleen (III_s) or both (IIIe+s)
- Stage IV Disseminated (multifocal) involvement of one or more extralymphatic organs, with or without associated lymph node involvement; or isolated extralymphatic organ involvement with distant (non-regional) nodal involvement

A and B Classification

Each stage should be divided into A and B according to the absence or presence of defined general symptoms. These are:

- 1. Unexplained weight loss of more than 10% of the usual body weight for the 6 months prior to first attendance.
- 2. Unexplained fever with temperature above 38 degree
- 3. Night sweat

Pathological Staging

Splenectomy, liver biopsy, lymph node biopsy and marrow biopsy are mandatory for the establishment of pathological staging

TNM classification of malignant tumours, Fifth Edition (1997)