Phone: 415-863-2380

Support: <u>info@pantomics.com</u>
Web: www.pantomics.com

Advancing Biomedical Science Through Tissue Amaga

Cat No: BRC15012 - Breast cancer tissue array

Lot#	Cores	Size	Cut	Format	QA/QC	1 2 3 4 5 6 7 8 9 10 11 12 13 14 15
BRC1501201	150	1.1mm	4um	10X15	H&E, IHC anti- Cytokeratin	Pantomics, In

Recommended applications: For Research use only. RNA or protein breast cancer/non-tumor tissue profiling using IHC or ISH; Antibody characterization.

Description: Breast cancer tissue array, set 12, non-overlapping with BRC1501 to BRC15011, 150 cores from normal-benign (5 cases) and cancer (70 cases) tissues in duplicates with grading, TNM staging data, AR, ER, PR, HER2, P53, EGFR, Ki67 IHC data.

All the tissues were from surgical resection. They were fixed in 10% neutral buffered formalin for 24 hours and processed using identical SOPs. Sections were picked onto Superfrost Plus or APES coated Superfrost slides. They can be stored for use at 4C for up to six months from the date of shipment. There may be 5 to 10% of tissue core loss.

Array Position	Sex	Age	Anatomic Site	Pathology	Grade	Stage
A01	F	52	Breast	Normal/hyperplasia	null	null
A02	F	49	Breast	Normal/hyperplasia	null	null
A03	F	39	Breast	Normal/hyperplasia	null	null
A04	F	22	Breast	Fibroadenoma	null	null
A05	F	23	Breast	Fibroadenoma	null	null
A06	F	50	Breast	Invasive ductal carcinoma	II	T2N0M0
A07	F	33	Breast	Invasive ductal carcinoma	II	T2N1M0
A08	F	56	Breast	Invasive ductal carcinoma	II	T2N0M0
A09	F	59	Breast	Invasive ductal carcinoma	II	T2N1M0
A10	F	37	Breast	Invasive ductal carcinoma	II	T2N0M0
A11	F	44	Breast	Ductal carcinoma in situ	null	TisN0M0
A12	F	46	Breast	Invasive ductal carcinoma	II	T3N0M0
A13	F	49	Breast	Invasive lobular carcinoma	II	T2N1M0
A14	F	42	Breast	Invasive ductal carcinoma	II	T3N1M0
A15	F	42	Breast	Invasive ductal carcinoma	II	T2N0M0
B01	F	52	Breast	Normal/hyperplasia	null	null
B02	F	49	Breast	Normal/hyperplasia	null	null
B03	F	39	Breast	Normal/hyperplasia	null	null
B04	F	22	Breast	Fibroadenoma	null	null
B05	F	23	Breast	Fibroadenoma	null	null
B06	F	50	Breast	Invasive ductal carcinoma	II	T2N0M0
B07	F	33	Breast	Invasive ductal carcinoma	II	T2N1M0
B08	F	56	Breast	Invasive ductal carcinoma	II	T2N0M0
B09	F	59	Breast	Invasive ductal carcinoma	II	T2N1M0
B10	F	37	Breast	Invasive ductal carcinoma	II	T2N0M0
B11	F	44	Breast	Ductal carcinoma in situ	null	TisN0M0
B12	F	46	Breast	Invasive ductal carcinoma	II	T3N0M0
B13	F	49	Breast	Invasive lobular carcinoma	II	T2N1M0
B14	F	42	Breast	Invasive ductal carcinoma	II	T3N1M0
B15	F	42	Breast	Invasive ductal carcinoma	II	T2N0M0
C01	F	34	Breast	Invasive ductal carcinoma	II	T2N0M0
C02	F	32	Breast	Invasive ductal carcinoma	III	T3N0M0
C03	F	41	Breast	Invasive ductal carcinoma	II	T2N0M0
C04	F	53	Breast	Invasive ductal carcinoma	II	T2N0M0
C05	F	44	Breast	Invasive ductal carcinoma	III	T2N0M0
C06	F	30	Breast	Invasive ductal carcinoma	II	T2N0M0
C07	F	36	Breast	Invasive ductal carcinoma	II	T2N0M0
C08	F	38	Breast	Invasive ductal carcinoma	II	T3N1M0
C09	F	28	Breast	Invasive ductal carcinoma	II	T2N0M0

Pantomics, Inc., 1049 Union Ave, Suite D, CA 94533 Phone: 415-863-2380; Fax: 510-225-1084; Email: info@pantomics.com

Phone: 415-863-2380 Support: <u>info@pantomics.com</u>

Web: www.pantomics.com

Pantomics, Inc.	anto		•		ww.panto	mics.com
		Advanci	ng Biomedical 9	Science Through Tissue Array		
C10	F	29	Breast	Invasive ductal carcinoma	II	T2N1M0
C11	F	54	Breast	Invasive ductal carcinoma	II	T2N0M0
C12	F	63	Breast	Invasive ductal carcinoma	II	T2N1M0
C13	F	56	Breast	Invasive ductal carcinoma	II	T2N0M0
C14	F	37	Breast	Invasive ductal carcinoma	II	T2N1M0
C15	F	76	Breast	Invasive ductal carcinoma	II	T2N0M0
D01	F	34	Breast	Invasive ductal carcinoma	II	T2N0M0
D01	F	32				
			Breast	Invasive ductal carcinoma	III	T3N0M0
D03	F	41	Breast	Invasive ductal carcinoma	II	T2N0M0
D04	F	53	Breast	Invasive ductal carcinoma	II	T2N0M0
D05	F	44	Breast	Invasive ductal carcinoma	III	T2N0M0
D06	F	30	Breast	Invasive ductal carcinoma	II	T2N0M0
D07	F	36	Breast	Invasive ductal carcinoma	II	T2N0M0
D08	F	38	Breast	Invasive ductal carcinoma	II	T3N1M0
D09	F	28	Breast	Invasive ductal carcinoma	II	T2N0M0
D10	F	29	Breast	Invasive ductal carcinoma	II	T2N1M0
D11	F	54	Breast	Invasive ductal carcinoma	II	T2N0M0
D12	F	63	Breast	Invasive ductal carcinoma	II	T2N1M0
D13	F	56	Breast	Invasive ductal carcinoma	II	T2N0M0
D14	F	37	Breast	Invasive ductal carcinoma	II	T2N1M0
D15	F	76	Breast	Invasive ductal carcinoma	II	T2N1M0
	F				II	T2N0M0
E01		43	Breast	Invasive ductal carcinoma		
E02	F	38	Breast	Invasive ductal carcinoma	III	T2N0M0
E03	F	50	Breast	Invasive ductal carcinoma	II	T2N1M0
E04	F	40	Breast	Invasive ductal carcinoma	II	T2N1M0
E05	F	47	Breast	Invasive ductal carcinoma	II	T2N0M0
E06	F	60	Breast	Invasive ductal carcinoma	II	T2N0M0
E07	F	56	Breast	Invasive ductal carcinoma	II	T2N0M0
E08	F	61	Breast	Invasive ductal carcinoma	II	T2N1M0
E09	F	38	Breast	Invasive ductal carcinoma	III	T2N1M0
E10	F	65	Breast	Invasive ductal carcinoma	II	T2N0M0
E11	F	61	Breast	Invasive lobular carcinoma	II	T2N0M0
E12	F	42	Breast	Invasive ductal carcinoma	II	T2N0M0
E13	F	24	Breast	Invasive ductal carcinoma	II	T2N0M0
E14	F	62	Breast	Invasive ductal carcinoma	II	T2N0M0
E15	F	54	+	Invasive ductal carcinoma	II	T2N0M0
	F	43	Breast		II	
F01			Breast	Invasive ductal carcinoma		T2N0M0
F02	F	38	Breast	Invasive ductal carcinoma	III	T2N0M0
F03	F	50	Breast	Invasive ductal carcinoma	II	T2N1M0
F04	F	40	Breast	Invasive ductal carcinoma	II	T2N1M0
F05	F	47	Breast	Invasive ductal carcinoma	II	T2N0M0
F06	F	60	Breast	Invasive ductal carcinoma	II	T2N0M0
F07	F	56	Breast	Invasive ductal carcinoma	II	T2N0M0
F08	F	61	Breast	Invasive ductal carcinoma	II	T2N1M0
F09	F	38	Breast	Invasive ductal carcinoma	III	T2N1M0
F10	F	65	Breast	Invasive ductal carcinoma	II	T2N0M0
F11	F	61	Breast	Invasive lobular carcinoma	II	T2N0M0
F12	F	42	Breast	Invasive ductal carcinoma	II	T2N0M0
F13	F	24	Breast	Invasive ductal carcinoma	II	T2N0M0
F14	F	62	Breast	Invasive ductal carcinoma	II	T2N0M0
F14 F15		54	+			
	F		Breast	Invasive ductal carcinoma	II	T2N0M0
G01	F	36	Breast	Invasive ductal carcinoma	II	T2N0M0
G02	F	38	Breast	Invasive ductal carcinoma	II	T2N0M0
G03	F	47	Breast	Invasive ductal carcinoma	III	T2N1M0
G04	F	46	Breast	Invasive ductal carcinoma	II	T2N0M0
G05	F	48	Breast	Invasive ductal carcinoma	III	T2N0M0
	F	59	Breast	Medullary carcinoma	null	T2N1M0
G06	!!					
G06 G07	F	60	Breast	Invasive ductal carcinoma	II	T2N0M0
		_	Breast Breast	Invasive ductal carcinoma Invasive ductal carcinoma	III	T2N0M0 T2N0M0
G07	F	60				

Pantomics, Inc., 1049 Union Ave, Suite D, CA 94533 Phone: 415-863-2380; Fax: 510-225-1084; Email: info@pantomics.com

Phone: 415-863-2380 Support: <u>info@pantomics.com</u>

Web: www.pantomics.com

			ng Biomedical S	cience Through Tissue Amay		
G11	F	36	Breast	Invasive ductal carcinoma	II	T2N0M0
G12	F	68	Breast	Invasive ductal carcinoma	II	T3N0M0
G13	F	49	Breast	Invasive ductal carcinoma	III	T2N1M0
G14	F	76	Breast	Invasive ductal carcinoma	II	T2N0M0
G15	F	77	Breast	Invasive ductal carcinoma	II	T2N0M0
H01	F	36	Breast	Invasive ductal carcinoma	II	T2N0M0
H02	F	38	Breast	Invasive ductal carcinoma	II	T2N0M0
H03	F	47	Breast	Invasive ductal carcinoma	III	T2N1M0
H04	F	46	Breast	Invasive ductal carcinoma	II	T2N0M0
H05	F	48	Breast	Invasive ductal carcinoma	III	T2N0M0
H06	F	59	Breast	Medullary carcinoma	null	T2N1M0
H07	F	60	Breast	Invasive ductal carcinoma	II	T2N0M0
H08	F	55	Breast	Invasive ductal carcinoma	III	T2N0M0
H09	F	58	Breast	Invasive ductal carcinoma	I	T2N0M0
H10	F	42	Breast	Invasive ductal carcinoma	II	T2N1M0
H11	F	36	Breast	Invasive ductal carcinoma	II	T2N1M0
H12	F	68	Breast	Invasive ductal carcinoma	II	T3N0M0
H13	F	49	Breast	Invasive ductal carcinoma	III	T2N1M0
H14	F	76	Breast	Invasive ductal carcinoma	II	T2N1M0
H15	F	77	Breast	Invasive ductal carcinoma	II	T2N0M0
I01	F	83	Breast	Invasive ductal carcinoma	II	T4N0M0
I02	F	46	Breast	Invasive ductal carcinoma	I	T2N0M0
I03	F	80	Breast	Invasive ductal carcinoma	II	T2N0M0
I04	F	43	Breast	Invasive ductal carcinoma	III	T3N0M0
I05	F	54			II	T2N3M0
106	F	59	Breast Breast	Invasive ductal carcinoma	II	T3N0M0
107	F	63		Invasive ductal carcinoma	null	
I08	F	51	Breast	Invasive papillary adenocarcinoma	III	T2N0M0 T2N0M0
108	F	55	Breast	Invasive ductal carcinoma	III	T2N0M0
		51	Breast	Invasive ductal carcinoma	III	
I10	F	58	Breast	Invasive ductal carcinoma	II	T2N0M0
I11	F	41	Breast	Invasive ductal carcinoma	I	T2N2M0
I12			Breast	Invasive ductal carcinoma		T1N1M0
I13	F	42	Breast	Invasive ductal carcinoma	II	T2N0M0
I14		72	Breast	Invasive ductal carcinoma	III	T2N0M0
I15	F	48	Breast	Invasive ductal carcinoma	II	T3N0M0
J01	F	83	Breast	Invasive ductal carcinoma	II	T4N0M0
J02	F	46	Breast	Invasive ductal carcinoma	I	T2N0M0
J03	F	80	Breast	Invasive ductal carcinoma	II	T2N0M0
J04	F	43	Breast	Invasive ductal carcinoma	III	T3N0M0
J05	F	54	Breast	Invasive ductal carcinoma	II	T2N3M0
J06	F	59	Breast	Invasive ductal carcinoma	II	T3N0M0
J07	F	63	Breast	Invasive papillary adenocarcinoma	null	T2N0M0
J08	F	51	Breast	Invasive ductal carcinoma	III	T2N0M0
J09	F	55	Breast	Invasive ductal carcinoma	III	T2N0M0
J10	F	51	Breast	Invasive ductal carcinoma	III	T2N0M0
J11	F	58	Breast	Invasive ductal carcinoma	II	T2N2M0
J12	F	41	Breast	Invasive ductal carcinoma	I	T1N1M0
J13	F	42	Breast	Invasive ductal carcinoma	II	T2N0M0
J14	F	72	Breast	Invasive ductal carcinoma	III	T2N0M0
J15	F	48	Breast	Invasive ductal carcinoma	II	T3N0M0

Notes: Bake at 60C for ~60 minutes before use. If antigen retrieving is needed, it is important to avoid **direct-boiling and high pH or high strength** antigen retrieving buffer. For availability of complimentary IHC data, please contact us at info@pantomics.com.

Phone: 415-863-2380

Support: <u>info@pantomics.com</u> Web: www.pantomics.com

Advancing Biomedical Science Through Tissue Arta

Certified by: Fancai Li, M.D.

TNM Classification: Breast carcinoma

T- Primary tumor

- TX Primary tumor cannot be assessed
- TO No evidence of primary tumor
- Tis Carcinoma in situ, intraductal carcinoma or lobular carcinoma in situ, or Paget's disease of the nipple with no tumor
- T1 Tumor 2 cm or less in greatest dimension
- T2 Tumor more than 2 cm but not more than 5 cm in greatest

dimension T3 - Tumor more than 5 cm in greasiest dimension

T4 - Tumor of any size with direct extension to chest wall or skin

N - Regional lymph nodes

- NX Regional lymph nodes cannot be assessed
- NO No regional lymph node metastasis
- N1 Metastasis to movable ipsilateral lymph node(s)
- N2 Metastasis to ipsilateral lymph nodes(s) fixed to one another or to other structure
- N3 Metastasis to ipsilateral internal mammary lymph node(s)

M - Distant metastasis

- MX Distant metastasis cannot be assessed
- MO No distant metastasis
- M1 Distant metastasis

TNM classification of malignant tumours, Fifth Edition (1997)